

Model 18007



WanderGuard[®] Options

Digital Timer

Model 18007

The digital timer enables facilities to program customized daily and weekly modes of operation for WanderGuard departure alert systems, WanderGuard Locks and Catch-All™ non-selective door alarms. The lighted LCD and keypad allow easy programming and viewing. Tamper resistant codes prevent unauthorized access or changes.



Model 18070

Passive Infrared (PIR) Monitor

Models 18060 and 18070

Passive Infrared (PIR) monitors allow expanded utilization of your WanderGuard departure alert system. This fully supervised device is easy to install for monitoring doors, hallways and elevators. Its advanced technology minimizes false alarms, while providing added safety for residents and staff.

Two models are available to meet your specific needs:

- Model 18060 for door and elevator applications
- Model 18070 for hallway applications

Model 18125



Bumper Guard – for WanderGuard monitors

Models 18125, 18140, and 18145

Protect your investment with WanderGuard bumper guards. These durable bumper guards are designed to protect your WanderGuard monitors from accidents and disasters. Made from **100% clear acrylic**, with a solid 1 inch thick base and 3/8 inch thick side panels, these cases blend into any decor.

- Model 18125 for WanderGuard E Series (8" x 3.5" x 2.75")
- Model 18140 for WanderGuard D Series (7.5" x 5.25" x 2.5")
- Model 18145 for FirstQ[®] Series (10.5" x 2.75" x 3.5")

STANLEY[®]

Healthcare Solutions

Toll-Free: 800-824-4478

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com