



Standard Features

- ▲ 4-door display
- ▲ 3 LEDs per monitored door
 - Green: Power indicator
 - Yellow: Pre-Alert alarm
 - Red: Departure alarm
- ▲ 2 distinct audible alarms for Pre-Alert and Departure
- ▲ Compatible with WanderGuard monitors and accessories

Specifications

- ▲ No separate power source required – powered by the door monitor (negative common)
- ▲ Dimensions: 10" x 2.5" x 1.5"
- ▲ Note: Model 14021 - The system installer should verify the included wiring meets local fire and building codes

Wired Remote Annunciator Systems

Wired 4-Channel Pre-Alert Remote Annunciator - Model 14021

The WanderGuard[®] 4-channel Remote Annunciator notifies staff when a resident wearing a WanderGuard signaling device is near, or exits through a WanderGuard-monitored door.

The patented Pre-Alert™ alarm display helps staff monitor elopements BEFORE they occur with a "condition yellow" LED display and alarm. If the resident exits, the Pre-Alert switches to a full departure alarm and a "condition red display".

Integration

The units can also be integrated with WanderGuard Locks and Catch-All™ non-selective door monitors; providing flexible, customized event notification at remote door locations not easily visible or assessable.

Nationwide Field Service

Our fleet of Field Support Representatives are company employees providing on-site installation, training and service.

STANLEY[®]

Healthcare Solutions

Toll-Free: 800-824-4478

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com



Standard Features

- ▲ 4-door display
- ▲ Single LED display per monitored door
- ▲ Audible alarm for departure notification
- ▲ Compatible with WanderGuard monitors and accessories

Specifications

- ▲ No separate power source required – powered by the door monitor (negative common)
- ▲ Dimensions: 10" x 2.5" x 1.5"
- ▲ Note: Model 14011 - The system installer should check for proper wiring which meets local fire and building codes

Wired Remote Annunciator Systems

Wired 4-Channel Remote Annunciator - Model 14011

The WanderGuard[®] 4-channel Remote Annunciator notifies staff when a resident wearing a WanderGuard signaling device exits through a WanderGuard-monitored door. The simple alarm display helps staff monitor elopements with a flashing LED and audible alarm.

Integration

The units can also be integrated with WanderGuard Locks and Catch-All™ non-selective door monitors; providing flexible, customized event notification at remote door locations not easily visible or assessable.

Nationwide Field Service

Our fleet of Field Support Representatives are company employees providing on-site installation, training and service.

STANLEY[®]

Healthcare Solutions

Toll-Free: 800-824-4478

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com