



Wireless Remote Annunciator Systems with Pre-Alert™ Model 14025

The **wireless 8-channel remote annunciator** offers the latest in remote event detection with digital 900 MHz electronics and **Pre-Alert** detection. Up to eight doors can be programmed with Pre-Alert and Alarm notification, each with distinct LEDs and audible alarm tones. **The patented Pre-Alert helps monitor egresses BEFORE they occur.** Supervisory indicators also notify you when the system is not functioning properly due to low battery, tampering or a missing transmitter. A numeric display indicates which transmitters and doors are being activated.

Standard Features

- ▲ **8-door programming capacity with Pre-Alert and Alarm**
- ▲ **Distinct audible alarms for Pre-Alert and Alarm notification**
- ▲ **Specific Pre-Alert and Alarm LEDs for each monitored door**
- ▲ **Three supervisory system LEDs for Tamper, Low Transmitter Battery and Inactive Transmitter**
- ▲ **Digital transmitter display**
- ▲ **Customized programming**
- ▲ **Large panel display**
- ▲ **Easy installation**
- ▲ **Writable label surface**
- ▲ **12 VDC plug-in power supply**

Interfacing Options

- ▲ **WanderGuard Departure Alert Systems**
- ▲ **WanderGuard Locks**
- ▲ **Catch-All non-selective door alarms**
- ▲ **TABS Select and Professional mobility monitors**
- ▲ **Arial wall transmitters**

Specifications

- ▲ **Power: Standard 120 VAC outlet**
- ▲ **Frequency: 900 MHz**
- ▲ **Dimensions: 11" x 8.5" x 1.75"**

Customized Programming

Allows for easy and custom programming to meet the specific needs of your facility. Transmitters are programmed directly from the base unit, and the display panel can be reprogrammed at any time.

Flexible Integration

The wireless 8-channel remote annunciator easily integrates WanderGuard® Locks, Catch-All™ non-selective door alarms, TABS® mobility monitors and Arial® wall transmitters with WanderGuard departure alert systems in the same panel to notify staff when a resident needs assistance.

Nationwide Field Service

Our fleet of Field Support Representatives are company employees providing on-site installation, training and service.

STANLEY[®]

Healthcare Solutions

Toll-Free: 800-824-4478

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com



Wireless Remote Annunciator Systems Model 14080

The 4-channel wireless remote annunciator offers an economical way to monitor a minimal number of event detectors from a remote location.

Flexible Integration

The wireless 4-channel remote annunciator easily integrates WanderGuard[®] Locks, Catch-All[™] non-selective door alarms, TABS[®] mobility monitors and Arial[®] wall transmitters with WanderGuard[®] departure alert systems in the same panel to notify staff when a resident needs assistance.

Standard Features

- ▲ Audible alarm
- ▲ Red light indicates alarm location
- ▲ Visual low battery indicator
- ▲ Visual power indicator
- ▲ Fully supervised (self-testing)
- ▲ Large panel display for labeling and identification
- ▲ 12 VDC plug-in power supply

Interfacing Options

- ▲ WanderGuard Departure Alert Systems
- ▲ WanderGuard Locks
- ▲ Catch-All non-selective door alarms
- ▲ TABS Select and Professional mobility monitors
- ▲ Arial wall transmitters

Specifications

- ▲ Power: Standard 120 VAC outlet
- ▲ Frequency: 900 MHz
- ▲ Dimensions: 10" x 2.5" x 1.5"

STANLEY[®]

Healthcare Solutions

Toll-Free: 800-824-4478

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com