


WanderGuard® Easy Programmer Model 18021

WanderGuard®
DEPARTURE ALERT SYSTEM



Overview

The Easy Programmer is a simple and effective way to program the WanderGuard ID signaling device bracelets with their own unique three digit identification number. The Easy Programmer will program the WanderGuard 90-day signaling device bracelet Model 36010. The WanderGuard ID Bracelet Manager is an on-board management database which stores all the bracelet data on both the Easy Programmer and the computer providing an automatic backup feature.

Features:

Easy Programmer includes:

- Unit, CD software, and USB connecting cable
- Assigns unique 3 digit identification numbers to each signaling device
- Reporting capabilities including signaling device expiration
- Integrates with Arial® Wireless Communication System

Specifications:

- Dimensions: 4" x 2.75" x 1.50"
- Weight: Approximately 3 oz.
- Power: Computer Enabled Power
- Computer Software Required: Windows XP Operating System with Service Pack 2 or higher
- USB connection: Required

STANLEY

Healthcare Solutions

Toll-Free: 800-824-4494

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com