

# TABS®

## Mobility Monitor



### Standard Features

- ▲ Limited Lifetime Warranty\*
- ▲ String and/or pressure pad monitoring
- ▲ Four distinct alarm tones
- ▲ Resets alarm during pressure pad use when patient returns his/her weight to the pad
- ▲ Low-battery alert
- ▲ Two wall mounting choices
- ▲ Locking garment clip
- ▲ No Velcro®

### Options

- ▲ Nurse call interface
- ▲ Bed bracket
- ▲ See Accessories page for more options

### Specifications

- ▲ Weight: Not more than 2 lbs
- ▲ Sound pressure level: >75dBA <90dBA
- ▲ Low battery indication
- ▲ Power Source: 9V alkaline battery
- ▲ Shunt Pull Force: Min. 10 oz; Max. 3 lbs

## TABS Select Monitor Model 25011

The TABS Select monitor is a restraint-free medical device that helps you:

- ▲ **Choose the appropriate monitoring methods for each individual: use the pull-string assembly, pressure pad or both simultaneously.**
- ▲ **Know when a patient is leaving a bed, chair, wheelchair, toilet or room.**
- ▲ **Identify a single resident or an entire wing with separate tones. Four distinct alarm tones are available.**
- ▲ **Avoid disturbances. When you add the optional nurse call interface, you can set the alarm to sound only at the nurse station and remain silent in the room.**

From the rugged case and secure mounting bracket to the locking garment clip, TABS products are some of the highest quality bed and chair exit alarms available.

Their versatility allows you to easily transfer a single monitor between bed, chair, wheelchair, and toilet.

Powered by a 9-volt battery, TABS monitors sound an alert to let you know when the batteries are low.

TABS monitors comply with JCAHO Patient Safety Standards. There are no on/off, volume control or delay switches, which helps to minimize human error.



locking garment clip



**Healthcare Solutions**

Toll-Free: 800-824-4490

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com