

Sidekicks ESA Bed Pad Model 26130

Provides the Extended Sensing Area needed for restless residents



FEATURES

- ▲ Restraint-free
- ▲ Maintains dignity and privacy
- ▲ Senses movement to help prevent falls
- ▲ 90-day limited warranty

SPECIFICATIONS

- ▲ Compatible with:
 - TABS Select (Model 25011)
 - TABS Professional (25022)
 - TABS Professional with Reset (25222)
- ▲ Dimensions:
 - 30" x 10" x .2"

NOTE: TABS Sidekicks pressure pads and other personal use products are not returnable.

Two Great Ways to Reduce the Risk of Falls

Our TABS® Professional and Select monitors act as the control unit for the TABS Sidekicks pressure pad system. All the great features of TABS Select and Professional monitors, such as multiple alarm tones and voice message alerts, can be used with the pressure pad.

This unique system gives you the flexibility to select the monitoring method that works best for each of your residents. Choose between the easy cord and clip assembly, or discreet pressure pads – or use them together – when you need to know when the resident is leaving a bed.

The bed model is placed on top of a mattress. The pad is then connected to the TABS monitor and the resident is situated comfortably on top of the pad. When the resident moves from the pad, the TABS monitor alarm sounds either with a local alarm or at the nurse call, or both.

TABS Sidekicks pressure pads and industry leading TABS mobility monitors combine to make the most versatile and cost effective fall monitoring system available.

