

TABS®

Mobility Monitor



Standard Features

- ▲ Restraint-free fall monitoring
- ▲ Works with the Arial wireless nurse call system
- ▲ Paging of fall-risk alerts
- ▲ Messages can be sent to multiple or individual pagers
- ▲ Resident identification
- ▲ Multiple alarm options
- ▲ Arial event tracking software

Arial® TABS Wireless Transmitter

Model 27000

A restraint-free mobility monitor for facilities equipped with the Arial wireless communication system

The Arial TABS wireless monitoring system integrates the benefits of the Arial nurse call system with TABS mobility monitors. This unique combination enables staff to respond to fall monitor alarms when they are not within hearing distance.

When an audible alarm is triggered on the TABS unit, a wireless signal is simultaneously sent to the Arial computer. The system also offers a paging option that allows staff to receive the TABS alarm message remotely within the facility while away from the nurse's station.

TABS monitors offer standardized training, easy-to-learn operation, and flexibility of use. Using a cord and clip assembly and/or TABS Sidekicks pressure pads, a single unit can be used on a bed, chair, or wheelchair. A special door mounting kit is also available enabling you to use a TABS device to alert staff when a patient opens the door to leave their room.

Arial TABS helps you create a comprehensive safety program by enhancing your Arial wireless nurse call system with effective fall management tools.