



Arial® SerWatch™ – Model 56050

Alerts staff to call system inactivity

The Arial wireless communication system utilizes the SerWatch feature to page staff if the system loses power, if the software or computer is interrupted, or if it is shut down. This “watchdog” feature gives facility managers additional assurance that staff members will receive resident calls. SerWatch monitors serial port activity and sends an “Arial Inactive” message to designated pagers if this activity is interrupted for any reason.

Standard Features

- ▲ Easily programmable through Arial software
- ▲ Monitors call system inactivity
- ▲ Notification of system interruption to staff pagers immediately

Benefits

- ▲ Enhanced paging reliability
- ▲ Immediate notification of computer or software interruptions
- ▲ Acts as a “watchdog” in the event of unexpected system down time

Specifications

- ▲ Software built into Arial paging base component

STANLEY

Healthcare Solutions

Toll-Free: 800-824-4494

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com