

Seat Belt Alarm Model 0707-945



FEATURES

- ▲ UMP monitor alarms when Velcro® closure is undone
- ▲ Adjustable strap offers a range of 18"-48"
- ▲ Fits most wheelchairs and side chairs
- ▲ Rugged construction
- ▲ Compatible with Standard Bed Sentry, Standard Chair Sentry, Deluxe Bed Sentry, Deluxe Chair Sentry, Voice Sentry, Deluxe Bed Alarm BS-5, Infrared Sentry IRM-1



A new way to reduce the risk of falls

The UMP® Seat Belt Alarm is designed to be used on wheelchairs to alert you if someone tries to exit the wheelchair. When the buckle is undone, the UMP monitor will alarm. The alarm message can be configured in multiple ways with all of the advanced features of the UMP Sentry™ line, such as multiple alarm tones or voice message alerts.

STANLEY

Healthcare Solutions

800-206-1045 | Fax 402-475-4023 | www.u-m-p.com