

Help Deter Resident Tampering

With New Dual Button Reset



Shown with optional cushion protector



Monitor Works with Revolutionary Timed Pads

Standard Chair Sentry™ Model 91610

Features:

- Helps deter resident tampering with:
 - Dual button reset; must be pressed simultaneously to silence alarm
 - Volume and delay switches inside battery compartment
 - No on/off switch
- Self-timing pad has visual and audible indicators for pad expiration
- Built-in Sentry Pad cord strain relief
- Low battery indicator
- Uses single 9V battery

Accessories:

- Complete line of Sentry Pads
 - Bed or chair
 - Self-timing
 - Moisture-resistant
 - Damage-resistant
 - Seat belt
 - 45-day, 90-day, 180-day, and 1-year options
- Latex-free cushion protector
- Universal bed bracket
- Universal chairstrap
- DC adapter

Distributed by:

Universal Medical Products

1.800.206.1045

1620 N 20th Circle

Lincoln, NE 68503

www.u-m-p.com

UIMP™

STANLEY

Revolutionary Timed Pads

Patent Pending



Features:

- Timer begins automatically after being plugged into equipped monitor
- LCD begins flashing when only 3 days remain and monitor emits a short audible alarm once per minute
- Number of days remaining in warranty period visually displayed
- Available in 45-day, 90-day, 180-day, and 1-year
- Available in bed and chair models

Full Line of Accessories:



Model 0440-137

Convenient T adapter for 1/4" mono plugs allows the use of push button call cord with monitor



Model 0707-569

Easily connect your monitor to your nurse call system with the nurse call cable. 8' with 1/4" mono jacks



Model 0707-566

Sturdy call button cord comes with clip for easy placement. 7 1/2' long with 1/4" mono jack



Model 0400-059

Keep your battery use down with this AC/DC adapter



Model 0707-578

Universal bed bracket works with all UMP metal brackets



Model 0707-752

Simple and easy to use the Universal Chairstrap allows all UMP monitors to be placed on oversized chairs