



From one door
to the whole floor

**we've got you
covered.**



**A fully scalable system
including call buttons, fall monitors,
and selective locking – all from one source.**

WanderGuard[®]
DEPARTURE ALERT SYSTEM



Trusted products & service for over 20 years.
Call 800-824-2996 or visit www.stanleyhealthcare.com

STANLEY[®]

Healthcare Solutions

A suite of products that can stand alone or integrate into a single comprehensive system.



WanderGuard® E Departure Alert System

Models 16516 and 16530

- Easily expanded to include resident identification, allowing for simple system upgrade as your facility grows
- Patented Pre-Alert™ detection helps monitor elopements before they occur
- Installs easily for use on single and double doors, hallways and elevators
- Compatible with WanderGuard D and 202 systems*
- Nationwide field service
- EMI alert
- Automatic self-test operation
- Power failure alert
- Digital keypad with programmable option
- Integrates with Arial® wireless communication system



WanderGuard ID Departure Alert System

Model 16901

- Resident ID on System Status Display
- Patented Pre-Alert detection helps monitor elopements before they occur
- Intelligent interference detection system
- Bracelet replace notification
- WanderGuard Bracelet Manager Software
- Installs easily for use on single and double doors, hallways and elevators
- Allows simple system upgrade and expandability as your facility grows
- Nationwide field service
- Automatic self-test operation
- Power failure alert
- Programmable keypad option
- Integrates with Arial wireless communication system



WanderGuard ID Departure Alert System

Model 16921

- Resident ID on System Status Display
- Patented Pre-Alert detection helps monitor elopements before they occur
- Intelligent interference detection system
- Bracelet replace notification
- WanderGuard Bracelet Manager Software
- Installs easily for use on single and double doors, hallways and elevators
- Allows simple system upgrade and expandability as your facility grows
- Nationwide field service
- Automatic self-test operation
- Power failure alert
- Programmable keypad option
- Integrates with Arial wireless communication system



Selective locking available with WanderGuard Locks

WanderGuard®
DEPARTURE ALERT SYSTEM

STANLEY®

Healthcare Solutions

800-824-2996 | Fax 402-475-4281 | www.stanleyhealthcare.com

*Only if these systems have been tested and/or returned for proper operation using WanderGuard E series signaling devices. ©2008-2010 Stanley Security Solutions, Inc. DOC-23-10964-AH