



Care starts with security

Give her the protection she deserves with the RoamAlert™ system

Reliable resident safety technology

STANLEY®

Healthcare Solutions

Elopement management solutions built for security

The RoamAlert system gives you the tools you need to protect your residents when and where they need it. No other product offers the same level of flexibility and scalability, combined with ease of use for staff and administrators.

Here's how the system works. Each resident wears a small radio transmitter (the smallest and lightest on the market). Exits are protected by door controllers.

When a resident approaches an exit, the door controller locks the door to prevent the resident from leaving; or if the door is open, an alarm sounds.

Beyond this essential functionality, the RoamAlert system offers you a wide range of options for central reporting, integration with other security systems, and even other applications such as personal emergency response and resident locating.

Security first

The RoamAlert system offers comprehensive security for individuals to walk about, based on proven technology.

The RoamAlert system's rapid response radio frequency technology means tags are *immediately* detected when at an exit, so that the door can be locked before the resident has a chance to open it.

A robust hardwired infrastructure ensures that system messages are always delivered to the server PC or other annunciation device. You never have to worry about messages disappearing into thin air.

And if the server PC is offline for some reason, your residents still remain protected: each door device is able to operate independently.

All system transactions are password controlled, time and date stamped and logged to a PC database. You have a complete record of exactly what happened, when and where, and which staff member responded.

Start small and grow

The RoamAlert system is modular and expandable. You can meet your needs today—and stay within your budget—while planning for the future.

Real-time resident locating, automatic door bypass, alarm reporting to portable devices... the RoamAlert system is capable of all this and more. But it can also be used to guard a single exit. And unlike other systems that require you to purchase an extensive infrastructure even for the most basic functionality, with RoamAlert you only have to buy what you need today.

When you're ready to expand, so is RoamAlert. The system can be scaled to cover not just an entire building, but an entire campus. You can add not just resident locating, but personal emergency response for your staff and residents. You can even track assets like the keys to medicine cabinets.



Built for your environment

The RoamAlert system offers a wealth of automatic features to ease the workload of caregivers, from monitoring alarms to escorting residents.

The RoamAlert system is simple to use, and staff require little training. All alarm information is presented in an intuitive visual format giving the resident's name, location and even a picture. Clearing alarms is quickly done at any PC or at the door.

For escorting residents, staff members simply need to enter their unique keypad passcode and then take the resident through the door. Or, you can assign each staff member a pendant: doors will unlock automatically.

The RoamAlert software includes several features that let you align system operation with your workflow. A Volunteer mode lets you assign certain resident tags to a pendant tag for bypass of selected exits by a volunteer. Finally, the system makes it easy to grant temporary absence to a monitored resident so they can pass through exits without generating an alarm. It's a secure and effective means to enable a family visit without sharing passcodes to bypass exits.

Resident centered

Every resident has unique needs. The RoamAlert system can help you design a program of care that is tailored to the individual.

There's no need for a 'one-size-fits-all' approach to elopement management. The RoamAlert system gives you flexibility to protect each resident as needed.

To allow the maximum freedom of movement, the system can be programmed to enable a particular resident to pass through certain exits—to reach a common area, say—while other exits remain protected.

For those residents who are particularly agitated and determined to remove a elopement management tag, the Securaband™ tag provides safe and comfortable protection. With a strap that resists cutting and a robust tag housing, the Securaband tag provides a tough physical barrier to unauthorized removal.



RoamAlert even helps you keep couples together, when one spouse begins to display wander symptoms. Their partner can be given a pendant that enables them to escort the resident—and that resident only—through any monitored exit.



Unique advantages of RoamAlert

False-alarm free technology

The RoamAlert system uses highly reliable FM technology and a hard-wired network that is not prone to interference from other radio frequency devices. No false alarms, and no missed alarms either.

Personalized resident protection

No other system offers as many options for customizing protection for each resident—letting them pass through some exits, but not others—giving them the maximum possible freedom.

Smallest resident tag

The RoamAlert tag is the smallest and lightest on the market today, only a third of an ounce and one inch across. With no hard edges, it is safe to use even on the most delicate skin, and so small that residents and family will barely notice it.

One system for all your security applications

The RoamAlert system can integrate a wide range of other systems and devices, from pagers to nurse call systems, giving you a single point from which to monitor activity around your facility.

Scalable and expandable

RoamAlert lets you start with one or two doors for perimeter protection and expand all the way to a facility-wide system for elopement management, emergency response and resident locating.

About Stanley Healthcare Solutions

Stanley Healthcare Solutions is a leading provider of safety, security, and supply chain management solutions for healthcare. Every day, thousands of hospitals and senior care facilities worldwide rely on Stanley products to deliver better care and improve efficiency.

Our smart storage solutions include supply/procedure carts, cabinets and casework, small carts, open storage, heavy-duty storage, and inventory management systems.

In the senior care market, Stanley Healthcare offers a full range of resident safety and security applications, including wander management, fall management and emergency call.

Our suite of patient security products provides individual protection to patients in multiple departments, from newborns in obstetrics to adults in the emergency department.

Stanley Healthcare Solutions will continue to expand its range of healthcare products and solutions for institutions and individuals, building on our years of experience in healthcare, and the Stanley legacy of quality and innovation.

For more information on the RoamAlert system and Stanley Healthcare Solutions' other products, please visit www.stanleyhealthcare.co.uk.