

Largest Selection of Pad Monitoring Options



Standard Pads

Features:

- Extra large size on bed pads helps to limit false alarms
- Durable vinyl cover cleans easily
- Compatible with all UMP® monitors
- Available in 45-day, 180-day, and 1-year



45-day Bed	92010
180-day Bed	92020
1-year Bed	92030
45-day Chair	93010
180-day Chair	93020
1-year Chair	93030



Damage Resistant Chair Pads

Features:

- Coiled cord design helps to eliminate tangled and pulled cords when used with a wheelchair
- Available in 90-day and 1-year

90-day Chair	93070
1-year Chair	93080



Moisture Resistant Chair and Bed Pads

Features:

- Innovative design eliminates the need for air holes
- Ultra thin design provides maximum comfort
- Compatible with all UMP monitors
- Available in 45-day and 90-day

45-day Bed	92050
90-day Bed	92060

45-day Chair	93090
90-day Chair	93100



Patent Pending

Timed Bed Pads

Features:

- Timer automatically begins after being plugged into equipped monitor
- Number of days remaining in usage period visually displayed
- LCD begins flashing when only 3 days remain and monitor emits a short audible alarm once per minute
- Available in 45-day, 180-day, and 1-year

45-day	92080
180-day	92100
1-year	92110



- Timer automatically begins when plugged into equipped monitor
- Timer starts at service life of pad and counts down to zero

**HELPS TO AVOID
F-TAGS**

Timed Chair Pads

Features:

- Timer automatically begins after being plugged into equipped monitor
- Number of days remaining in usage period visually displayed
- LCD begins flashing when only 3 days remain and monitor emits a short audible alarm once per minute
- Available in 45-day, 180-day, and 1-year



Patent Pending

45-day	93120
180-day	93140
1-year	93150