

Data Sheet:

## TABS® Sensor Pads

Two Great Ways to Reduce the Risk of Falls

### Overview

TABS Professional monitors act as the control unit for the TABS sensor pad system. All the great features of TABS Professional monitors, such as multiple alarm tones and voice message alerts, can be used with the sensor pad.

This unique system gives you the flexibility to select the monitoring method that works best for each of your residents. Choose between the easy cord and clip assembly, or discreet pressure pads—or use them together—when you need to know when the resident is leaving a bed, chair or wheelchair.

The chair model can be placed on a wheelchair or other seat. The bed model is placed on top of a mattress and underneath the fitted sheet. The pad is then connected to the TABS monitor and the individual is situated comfortably on top of the pad. When the individual moves from the pad, the TABS monitor alarm sounds either with a local alarm or at the nurse call or both.

TABS sensor pads and industry-leading TABS mobility monitors combine to make the most versatile and cost effective fall monitoring system available.

NOTE: TABS sensor pads and other personal use products are not returnable.

### Models 26000 and 26100



### Product Highlights

- Restraint-Free
- Maintains individual's dignity and privacy
- Senses movement to help reduce falls
- 90-day limited warranty

### Specifications

#### Control Units

- Tabs Professional Monitor
- Power Source: 9V alkaline battery

#### Dimensions

- Chair Model 26000: 7.0" x 15.0" x 0.15" (17.78 cm x 38.1 cm x 0.38 cm)
- Bed Model 26100: 10.0" x 30.0" x 0.15" (25.4 cm x 76.2 cm x 0.38 cm)

STANLEY Healthcare  
130 Turner Street, Waltham, MA 02453

Phone: +1-888-622-6992  
Email: stanleyhealthcare@sbdinc.com

www.stanleyhealthcare.com