

# TABS®

## Mobility Monitor



### Features

- ▲ Lifetime warranty\*
- ▲ Audible initial pressure detection
- ▲ Separate sensor pad input jack
- ▲ Separate nurse call input jack
- ▲ String and/or pressure pad monitoring
- ▲ Built-in microphone
- ▲ Four distinct alarm tones
- ▲ Resets alarm during pressure pad use when patient returns weight to pad
- ▲ Low-battery light and alert
- ▲ Two wall mounting choices
- ▲ Locking garment clip
- ▲ No Velcro®
- ▲ ETL® Listed
- ▲ Latex free

### Options

- ▲ Nurse call interface cord
- ▲ AC power supply
- ▲ Wire bracket
- ▲ See Accessories page for more options

### Specifications

- ▲ Weight: Not more than 2 lbs
- ▲ Sound Pressure level: >75dBA <90dBA
- ▲ Power Source: 9V Alkaline Battery
- ▲ Shunt Pull Force: Min. 10 oz; Max. 3 lbs

\* Contact Stanley Healthcare Solutions at 800-824-4490 for a copy of the TABS Limited Lifetime Warranty in its entirety.

## TABS Professional Monitor Model 25022

The TABS Professional monitor is a bed and chair exit monitor that helps you improve patient safety by allowing you to:

- ▲ Easily use nurse call and sensor pads simultaneously due to the separate input jacks.
- ▲ Remind patients to call for help or use their cane or walker.
- ▲ Easily record messages or customize alarm tones with the built-in microphone.
- ▲ Choose the appropriate monitoring method for each individual whether pressure pad, pull-string assembly or both.
- ▲ Avoid disturbing other patients. By adding the optional nurse call interface, you can set the alarm to sound only at the nurse station or have a voice message at the monitor and an alarm at the nurse station.

From the rugged case and secure mounting bracket to the locking garment clip, TABS products are some of the highest quality bed and chair exit monitors available.

Their versatility allows you to easily transfer a single monitor between bed, chair, wheelchair and toilet.

Powered by a 9-volt battery (optional AC power supply available), TABS monitors sound an alert with a light to let you know when the batteries are low.

The TABS Professional is **latex free** and is **ETL® listed to comply with UL® 1069**. TABS monitors also **support JCAHO Patient Safety Standards**. There are no on/off, volume control or delay switches, which helps to minimize human error.



**STANLEY®**  
**Healthcare Solutions**

Toll-Free: 800-824-4490

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com