

TABS®

Mobility Monitor



Standard Features

- ▲ Lifetime warranty*
- ▲ Easy to use and install
- ▲ Durable, ergonomic design
- ▲ Easy transfer between bed, wheelchair and room
- ▲ Two wall mounting choices
- ▲ 30-day return privilege
- ▲ Locking garment clip
- ▲ No Velcro®

Options

- ▲ Bed bracket
- ▲ Over-the-door bracket
- ▲ Various pull cord lengths

Specifications

- ▲ Weight: Not more than 2 lbs
- ▲ Sound Pressure level: >75dBA <90dBA
- ▲ Low Battery Indication
- ▲ Power Source: 9V Alkaline Battery
- ▲ Shunt Pull Force: Min. 10 oz; Max. 3 lbs

TABS Basic Monitor Model 25000

The basic tool for fall management

The TABS mobility monitor is a restraint-free medical device that:

- ▲ Alerts you when a resident or patient is leaving a bed, chair, wheelchair or room
- ▲ May serve as a reminder to residents or patients to call for assistance before rising
- ▲ Helps to reduce the risk of falls
- ▲ Helps improve staff response time

With a simple disk, clip, and cord system, the TABS Basic monitor allows you to determine a personal safety zone for each resident and then alerts you when residents leave their personal safety zones.

TABS Sidekicks pressure pads are not compatible with the TABS Basic model. There are other TABS monitors available that can use sensor pads. Call your TABS sales representative to determine which monitor is best suited for your purpose.

From the rugged case and secure mounting bracket to the locking garment clip, TABS products are among the highest quality fall prevention monitors available. Their versatility allows you to easily transfer a single monitor between bed, chair, wheelchair, toilet, and door. The over-the-door bracket enables TABS alarms to be used as a room monitor.

Powered by a 9-volt battery, TABS monitors sound an alert to let you know when the batteries are low.

TABS monitors comply with JCAHO Patient Safety Standards. There are no on/off, volume control or delay switches, which helps to minimize human error.



locking garment clip



Healthcare Solutions

Toll-Free: 800-824-4490

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com