

TABS® Accessories

Optional Cord Lengths



Optional Cord Lengths

TABS monitors are shipped with a 30 inch cord. The cord, however, must be long enough to allow motion within the caregiver's predetermined safety zone, but short enough to ensure that the alarm is activated when the patient or resident moves beyond the safety zone. Different placements of the TABS monitor may require different lengths of cord. Make sure you'll have the proper length of cord and replacement cords on hand when you need them by ordering in advance.

- | | | | |
|-----------------------|---|----------------------------------|---|
| Model 24040: 15" cord | ○ | Model 24050: 18" cord | ● |
| Model 24060: 24" cord | ● | Model 24070: 30" cord (standard) | ● |
| Model 24080: 36" cord | ● | Model 24090: 42" cord | ● |

TABS pull cords feature a locking garment clip.

Cordless Reset Option



Cordless Reset Option Model 0707-280

When using TABS monitors with pressure pads, you may not want to use the cord and clip assembly. The TABS cordless reset option is a metal shunt that replaces the cord and clip assembly to keep the TABS monitor from alarming without the additional cord.

Nurse Call Interface



Nurse Call Interface Model 25150

TABS Professional and Select monitors interface easily with most call systems with this simple adapter.



Healthcare Solutions

Toll-Free: 800-824-4490

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com

TABS® Mounting Options

Dual Lock® Fastener Mounting



Dual Lock Fastener Model 24130

On surfaces where more permanent mounting with screws is impractical, use Dual Lock fasteners to mount your TABS monitor. Unlike Velcro®, Dual Lock can be disinfected. Not recommended for wallpapers or veneers.

One piece measures 2"x 2 1/2"

Wall Mounting Brackets



Wall Mounting Bracket Model 25120

To get the most out of your TABS mobility monitor, install mounting brackets wherever the resident or patient will need monitoring. Add a bracket to the wall behind a favorite chair, behind the commode or even at a loved one's home.

Wheelchair Brackets



Easy-Clip Wheelchair Bracket Model 25300

The Easy-Clip bracket can be used on standard wheelchairs and Merry Walkers®. Installing in seconds, it enables the TABS monitor to effortlessly snap into place on the secure bracket.

Bed Bracket



Bed Bracket Model 0120-123

The TABS bed bracket enables you to temporarily attach a TABS monitor to most headboards. The secure bracket can be easily removed from the bed without damaging the headboard's veneer.

STANLEY®

Healthcare Solutions

Toll-Free: 800-824-4490

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com