

Seat Belt Alarm Model 24201



FEATURES

- ▲ TABS monitor alarms when Velcro® closure is undone
- ▲ Adjustable strap offers a range of 18"-48"
- ▲ Fits most wheelchairs and side chairs
- ▲ Rugged construction
- ▲ Compatible with:
 - TABS Select (Model 25011)
 - TABS Professional (25022)
 - TABS Professional with Reset (25222)



A new way to reduce the risk of falls

The TABS® Seat Belt Alarm is designed to be used on wheelchairs to alert you if someone tries to exit the wheelchair. When the buckle is undone, the TABS monitor will alarm. The alarm message can be configured in multiple ways with all of the advanced features of the TABS Select or TABS Professional, such as multiple alarm tones, voice message alerts or wireless integration into the Arial® communication system for wireless paging.

STANLEY

800-824-4490 | Fax 402-475-4281 | www.stanleyhealthcare.com

Healthcare Solutions