



Standard Features

- ▲ Fully supervised (self-testing)
- ▲ Low battery alert
- ▲ User-replaceable, long-life lithium battery
- ▲ Water-resistant pendant transmitter
- ▲ Braille lettering on wall transmitter button
- ▲ Adjustable length pullcord for wall transmitter
- ▲ Easy to install slide and lock wall bracket
- ▲ Red light flashes when a call is sent
- ▲ Reset / check-in button feature built in

Benefits

- ▲ Single or double push button cords for monitoring two residents (wall transmitter only)
- ▲ Wristband (pendant transmitter)
- ▲ Additional wall mounting brackets
- ▲ Wireless locator to identify approximate location of pendant/transmitter calls

Specifications

- ▲ Power: 3V lithium battery
- ▲ Frequency: 900 MHz spread spectrum, 100 KHz bandwidth
- ▲ Wall Transmitter
Dimensions: 3.4" x 5" x 1.4"
Weight: 4.8 oz.
- ▲ Pendant Transmitter
Dimensions: 1.6" x 2.9" x 7.2"
Weight: 1.4 oz.

Arial® Transmitters

Arial transmitters enable your residents to place calls to staff from anywhere within the system coverage area. Residents gain independence to pursue their own interests, and the peace of mind knowing help is only a push button away.

The Arial system offers versatility. The lightweight, water-resistant pendant transmitter can be worn as a necklace, on a belt, or with an optional wristband. Wall transmitters include slide and lock brackets. Installing additional brackets makes it easy for you to utilize a single transmitter in multiple locations.

With the Arial system, a call for assistance is activated by pressing the extra-large wall transmitter button, by pulling the adjustable length pull cord, or by pressing the pendant transmitter button. You may also use optional push-button cords with the wall transmitter for non-ambulatory residents.

The check-in button on the 54036 wall transmitter enables residents to send an "all-OK" message. This allows staff to quickly account for all residents with a daily report.

You'll feel more confident knowing Arial is a fully supervised system that tests itself and will alert you if a transmitter is not functioning or has a low battery. Arial is an expandable system that allows you to economically add transmitters as your facility grows.

Available Models:

Wall Transmitter with Check-in
Model 54036

Water Resistant Wall Transmitter
Model 54057

Water Resistant Pendant
Model 59050

Pendant Options:

Necklace clip (included)

Belt clip (included)

Wristband clip
8 inch - 0707-205

10 inch - 0707-220

Pendant Accessory Kit
0707-586
(Chain, clip, and belt clip)

STANLEY®

Healthcare Solutions

Toll-Free: 800-824-4494

Office: 402-475-4002

Fax: 402-475-4281

www.stanleyhealthcare.com