

UIMP High-Tech, Low Cost Stanley® Solutions

Optional Accessories



Shown with Optional
Cushion Protector
Model 0141-409



Universal Bed Bracket
Model 0707-578



Economy Pull String Model 91230

Features:

- Easy-to-use, low cost, pull-string monitor
- Single person monitoring device that alarms when unsupervised resident movement occurs
- Monitor attaches to a bed or chair; Safety Klip is attached to the resident
- Alarm sounds when magnetic seal is broken between magnet and monitor, indicating resident movement

Safety Klip

This patent-pending device is locked onto an article of the resident's clothing. The unique plastic tab makes it difficult for the resident to detach without the assistance of a caregiver.

Distributed by:

Universal Medical Products
1.800.206.1045
1620 N 20th Circle
Lincoln, NE 68503
www.u-m-p.com

