

# U-M-P High-Tech, Low Cost Stanley® Solutions

## Optional Accessories

### Cushion Protector



Model 0141-497

### Nurse Call Cable



Model 0707-569

### Universal Bed Bracket



Model 0707-578

### Call Cord



Model 0707-566

### T-connector for Nurse Call



Model 0440-137



## Deluxe Pull String Model 91650

### Features:

- Pull-string magnet monitor that alarms when unsupervised resident movement occurs
- Attaches to a bed or chair, and the exclusive "Safety Klip" is attached to the resident
- Alarm sounds when the magnetic seal is broken between the monitor and the magnet
- Can be integrated into nurse call system

### Safety Klip

This patent-pending device is locked onto an article of the resident's clothing. The unique plastic tab makes it difficult for the resident to detach without the assistance of a caregiver.

Distributed by:

Universal Medical Products  
**1.800.206.1045**  
1620 N 20<sup>th</sup> Circle  
Lincoln, NE 68503  
www.u-m-p.com

