

When safety is your first concern,
quality counts.

Tabs[®]
MOBILITY MONITOR



**Because every patient and resident is unique,
we offer versatile products that can adapt
to the needs of each individual.**

TABS[®] products allow you to create individualized care plans that will help you meet F-Tag 323 expectations on Adequate Supervision and use of safe equipment.

We'll design and customize a system that meets your requirements and budget with the quality only offered by the makers of WanderGuard[®].



From the
makers of
WanderGuard

Trusted products & service for over 20 years.
800-824-2996 www.stanleyhealthcare.com

STANLEY[®]

Healthcare Solutions

Quality Fall Management Tools to help you focus on individualized care.



TABS Basic Monitor

Model 25000

- Simple to use cord and clip monitoring
- Easy transfer between bed, wheelchair and room
- Low battery alert for peace of mind



TABS Select Monitor

Model 25011

- Patented dual or selective monitoring with cord and clip and TABS Sidekicks pressure pad
- Easy transfer between bed, wheelchair and room
- Four distinct alarm tones distinguish monitored patients
- Can be connected to your existing nurse call system
- Low battery alert for peace of mind



TABS Professional Monitor

Model 25022

- Audible initial pressure detection
- 110 Volt power capability
- String and/or pressure pad monitoring with TABS Sidekicks
- Recordable voice message with built-in microphone
- Separate sensor pad input jack
- Four distinct alarm tones
- Resets alarm during pressure pad use when patient returns their weight to the pad
- Two wall mounting choices
- Low battery alert for peace of mind
- Can be connected to your existing nurse call system



TABS Elite Monitor

Model 29800

- Discreet bed monitoring with TABS Sidekicks pressure pad
- Can be connected to most existing nurse call systems
- Volume control with silent alarm option (connected to nurse call)
- Time delay control for restless or sedate patients
- Powered by standard 120V outlet



TABS Sidekicks Pressure Pad

Models 26000 and 26100

- Restraint-free, discreet monitoring maintains dignity and privacy
- Senses movement to help prevent falls
- Use on chair or bed



Locking Garment Clip

Available on all models



TABS
MOBILITY MONITOR

STANLEY

800-824-2996 | Fax 402-475-4281 | www.stanleyhealthcare.com

Healthcare Solutions

*Contact us for a copy of the TABS Limited Lifetime Warranty in its entirety. ©2008-2010 Stanley Security Solutions, Inc. DOC-23-25034-AC