

...I want a customizable, integrable call system that's versatile yet easy to use. Where do I find that?



ARIAL[®]

Wireless Communication System

Arial helps you provide a higher caliber of care with flexible nurse call options specific to your needs and your budget, with the freedom to expand over time.

Arial nurse call devices constantly report their status to the Arial system. Staff members will be informed if a battery runs low, if a transmitter isn't working properly, or if a transmitter leaves the coverage area.



Trusted products & service for over 20 years.
Call 800-824-2996 or visit www.stanleyhealthcare.com

STANLEY[®]

Healthcare Solutions

Reducing risk, increasing safety and optimizing staff time needs more than an emergency call system.



Arial Wireless Communication System

- Nurse call system for personal monitoring needs of over 1,000 patients/residents
- Receives calls for help from wall transmitters, personal pendants, or blow tubes at the central monitoring station and at pagers carried by staff
- State-of-the-art 900MHz wireless technology
- Integrates with smoke detectors, door/window transmitters, passive infrared detectors, and other personal monitoring devices from Stanley Healthcare Solutions
- Records call history to help managers develop better protocols for quicker response times
- Can augment wired systems to accommodate your facility's growth
- Integrable with other Stanley-Senior Technologies products to achieve additional monitoring, e.g. wanderer management, fall monitoring, and locking



Arial Transmitters

Models 54036, 54037, 54057 and 59050

- Wall transmitters available with check-in or reset features built in
- Water resistant wall and pendant models available
- Fully supervised (self-testing)
- Low battery alert
- User-replaceable, long-life lithium battery
- Braille lettering on wall transmitter button
- Adjustable length pullcord for wall transmitter
- Easy to install slide and lock wall bracket
- Red light flashes when a call is sent
- Wireless locator to identify approximate location of pendant/transmitter calls
- Necklace and belt clips included with pendant; optional wristband clip and accessory kit available
- Expandable system allows additional transmitters as your facility grows



Arial Paging Options

Models 52111, 52112 and 52113

- The Apollo features a mid-sized case with a larger screen, covered LED and increased durability
- 14 selectable alert tones
- Unread message notification
- Silent vibrating alert
- Backlit display for low light environments
- 26-slot message capacity
- 16 month, date, year, time display
- Zoom feature for bigger text
- Daily alarm set
- Duplicate message detection
- Programmable for 2 or 4 line messages
- Total character memory: 32,000
- Uses one AAA battery
- Available with 30-day, 90-day, or 1-year warranty
- Paging zones can be customized specifically for your facility

Integrates with WanderGuard®, TABS® and Bed-Check® products



ARIAL®

Wireless Communication System

STANLEY®

Healthcare Solutions

800-824-2996 | Fax 402-475-4281 | www.stanleyhealthcare.com