

Accessories



Model 0707-752

Simple and easy to use, the universal Chair Strap allows all UMP monitors to be placed on oversized chairs.



Model 0400-059

Keep your battery use down with this AC/DC adapter.



Model 0440-137

Convenient T-adaptor for 1/4" mono plugs allows the use of push button call cord with monitor.



Model 0707-569

Easily connect your monitor to your nurse call system with the nurse call cable. 8' with 1/4" mono jacks.



Model 0707-566

Sturdy call button cord comes with clip for easy placement 7 1/2' long with 1/4" mono jack.



Model 0707-754

Receiver extender allows for monitor placement away from resident accessible areas such as behind the headboard.



Model 0707-578

Universal bed bracket works with all UMP fall monitors.



Model 0707-755

Resets alarm at a distance. Integrates with most name badges. Each Remote Reset Pendant works with all properly equipped UMP monitors. **Patent Pending**