

# Deluxe Bed Sentry



## Deluxe Bed Sentry

### Features:

- State-of-the-art digital electrical Bed Monitoring System with exclusive 2' x 3' Air Sheet; the Monitor mounts to the headboard with the Air Sheet inserted beneath the mattress
- Alarm is reset by returning the resident to bed or using the caregiver hold button
- Silent alarm option sounds alert only at nurse's station
- Works with 2-year Air Sheet
- 0 - 6 second alarm delay controls
- Nurse call interface
- Powered by A/C adapter (included)

Model	95360
-------	-------