

# UMP® Fall Management

## HELPS PROTECT YOUR RESIDENTS AND REDUCES FALL RISK

To provide high-quality care, senior communities nationwide are relying on STANLEY Healthcare's UMP Fall Management solutions to monitor residents who are at risk of falling, whether in a bed, a wheelchair, or the restroom.

To reduce the common challenge of alarm fatigue, UMP monitors can be integrated with your existing nurse call and emergency call systems, such as STANLEY Healthcare's Arial® Emergency Call solution.



### UMP SENTRY™ MONITOR MODEL 91621

- Dual button reset—buttons are pressed simultaneously to silence alarm
- Volume and delay switches inside battery compartment
- No on/off switch
- Built-in Sentry pad cord strain relief
- Low battery indicator—uses single 9V battery
- Selectable 4-tone alarm
- Optional nurse call integration
- AC/DC adapter jack
- Integrated cushion protector
- Hinged battery cover



### SENSOR PAD MONITORING OPTIONS

- Durable vinyl cover for easy cleaning
- No holes construction

**45-day Bed** Model 92010

**180-day Bed** Model 92020

**1-year Bed** Model 92030



**45-day Chair** 93010

**180-day Chair** 93020

**1-year Chair** 93030



### ACCESSORIES

<b>Universal Bed Bracket</b>	Model 0707-578
<b>AC/DC Adapter</b>	Model 0400-059
<b>Seat Belt</b>	Model 0707-945
<b>Floor Mat</b>	Model 93040
<b>Arial ES Transmitter</b>	Model 54351
<b>Nurse Call Button Cord</b>	Model 0707-566
<b>Nurse Call Cable</b>	Model 0707-569