

SOLUTION OVERVIEW

# M200 Fall Monitoring System

Supporting safety and individualized care



## Solution Benefits

- **Address many needs with one solution.** No need to purchase and maintain multiple solutions to meet the diverse needs of the individuals in your care. Choose the delay times, tamper resistance and alert type settings that best support each person's situation and care plan.
- **Give caregivers more time to focus on people.** The M200 reduces the need to continually disconnect and reconnect pad cords, and will help reduce nuisance alarms, alarm fatigue, and enable caregivers to spend more time on high-value interactions with patients or residents.
- **Simplify maintenance and operation.** Because it can be configured to meet specific needs, the M200 makes it possible for standardization to a single solution throughout your organization. That reduces training requirements and streamlines ongoing maintenance and operation.

## Challenges

Implementing and maintaining a fall management program to reduce the quantity and severity of falls. Developing and supporting individual care plans, including fall prevention measures.

## Solution

The M200 Fall Monitoring System is many solutions in one. You can configure it to support your unique fall management program—and to address the individual needs of each patient or resident.

## How It Works

With the M200 Fall Monitoring System, you can connect two pads (one on a bed, another on a wheelchair or side chair) to a single monitor. Select the delay setting (0, 1 or 2 seconds) that best protects each person in your care. Choose the configuration that best suits the criteria to successfully execute care plans and program goals: electric or battery power, auto-reset or non-auto-reset, audible or silent alerts, and tamper-resistant or easy-access rear cover.



*The M200 Fall Monitoring System integrates with virtually any nurse call or e-call solution—positioning it to become an indispensable part of your fall management program.*



### About STANLEY Healthcare

Over 15,000 hospitals and senior living communities rely on STANLEY Healthcare solutions to empower caregivers to deliver better care. STANLEY Healthcare is a part of Stanley Black & Decker and a proud supporter of the Alzheimer's Association®. Learn more at [stanleyhealthcare.com](http://stanleyhealthcare.com).