

Data Sheet:

Sentry™ Model 91621

Help Deter Resident Tampering

Features

- Helps deter resident tampering with:
 - Dual button reset; must be pressed simultaneously to silence alarm
 - Volume and delay switches inside battery compartment
 - No on/off switch
- Wireless remote reset option
- Self-timing pad has visual indicator for pad expiration
- Built-in Sentry Pad cord strain relief
- Low battery indicator
- Uses single 9V battery
- Selectable 4-tone alarm
- Nurse call with remote alarm option
- DC adapter jack

Optional Accessories

- Complete line of sensor pads available for bed or chair
- Floor mat available
- Integrated Latex-free cushion protector
- Universal bed bracket
- Universal chairstrap

