



COVID-19 QUICK DEPLOYMENT BUNDLE

Asset Management with ICU Medical Integration

Challenges

Due to the recent COVID-19 outbreak, it is essential for hospitals and nursing departments to quickly locate and identify the status of critical equipment such as infusion pumps. Understanding the status of a pump (e.g., In-use, Available) allows caregivers to quickly locate ones that are ready for use and to increase the utilization of these critical assets.

Solution

STANLEY Healthcare's AeroScout® Asset Management solution is a Real-Time Location System (RTLS) that automatically monitors and manages the location and status of critical equipment throughout hospitals and healthcare facilities.

By combining your ICU Medical MedNet™ IV safety software with STANLEY Healthcare's AeroScout® Real-Time Location System (RTLS) platform, you gain a bi-directional connection that provides real-time visibility into your IV pumps' location and status for clinical use and maintenance, helping you save time while enhancing overall fleet operation.

How It Works

An RTLS tag is securely attached to your ICU Medical Plum 360™ and LifeCare PCA™ infusion pumps, providing location information for the device. The integrated communication network between the STANLEY Healthcare AeroScout platform and ICU Medical MedNet™ Software provides real-time status information—*In-Use*, *Not-In-Use* or *Standby*—to help efficiently locate pumps for use, maintenance or distribution while bypassing pumps in clinical use. Access to this information is easily viewed in either system for added workflow flexibility and convenience.

Deployment

We have designed a simplified and remote implementation process that enables hospitals and healthcare facilities to quickly take advantage of our RTLS integration within the Asset Management solution. This enables you to be fully operational on an accelerated schedule.

COVID-19 RESPONSE LIMITED-TIME OFFER*


To support the needs of our customers as they work to reduce the spread of COVID-19 while caring for infected patients, we have developed a special package that enables you to quickly deploy the ICU Medical infusion pump Integration for locating and managing the location and status of critical infusion pumps.



The Value of Integration

Actively used by multiple hospitals, this bi-directional solution provides status and location information of the Plum 360 and LifeCare PCA 7 devices for faster, more efficient pump management.

The integrated solution provides valuable information to streamline biomedical engineering workflows and reduce staff search time and labor costs. Meanwhile, sterile processing professionals can quickly understand what pumps require cleaning and redistribution within a PAR-managed inventory process. And, clinicians can feel confident in their ability to quickly and easily find the nearest available infusion pump for clinical use.



Visit our **COVID-19 Resource Center** for updated information to help you manage the pandemic.

www.stanleyhealthcare.com/covid-19

CONTACT US

Contact your STANLEY Healthcare Client Executive to learn more, or call **+1-888-622-6992**

*DISCLAIMER

All product and pricing offerings are subject to a mutually agreed quote and acceptance by STANLEY Healthcare of customer's purchase order. STANLEY Healthcare reserves the right to set criteria for this promotion based on factors including, but not limited to, contractual terms and product availability. The offer period begins on April 1, 2020 and expires on June 30, 2020. Customer must submit a purchase order to STANLEY Healthcare no later than June 30, 2020. This offer may not be applied retroactively for orders placed prior to April 1, 2020. STANLEY Healthcare reserves the right to modify or cancel this offer at any time without notice. Offer valid for customers in the United States and Canada only.

STANLEY
Healthcare

About STANLEY Healthcare

Over 15,000 hospitals and senior living communities rely on STANLEY Healthcare solutions to empower caregivers to deliver better care. STANLEY Healthcare is a part of Stanley Black & Decker and a proud supporter of the Alzheimer's Association®. Learn more at stanleyhealthcare.com.

icumedical
human connections

About ICU Medical

ICU Medical connects patients and caregivers through safe, life saving, life enhancing IV therapy products, systems, and services. Since IV therapy is our only business, meeting your IV needs with quality products and consistent supply is our only concern. To learn more, visit icumed.com.

75 Portsmouth Blvd., Suite 220, Portsmouth, NH 03801 | 1-888-622-6992 | stanleyhealthcare.com | stanleyhealthcare@sbdinc.com