

SOLUTION OVERVIEW

M210 Fall Monitoring System

Supporting safety and individualized care



Solution Benefits

- **Address many needs with one solution.** No need to purchase and maintain multiple solutions to meet the diverse needs of the individuals in your care. Choose the delay times, tamper resistance and alert type settings that best support each person's situation and care plan.
- **Give caregivers more time to focus on people.** The M210 reduces the need to continually disconnect and reconnect pad cords, and helps reduce nuisance alarms, alarm fatigue, and enable caregivers to spend more time on high-value interactions with patients or residents.
- **Simplify maintenance and operation.** The M210 firmware can be updated with the new mobile app, eliminating the need to replace monitors every couple of years. Standardizing on a single solution reduces training requirements and streamlines ongoing maintenance and operation.

Challenges

Implementing and maintaining a fall management program to reduce the quantity and severity of falls. Developing and supporting individual care plans, including fall prevention measures.

Solution

The M210 Fall Monitoring system is an enhanced version of our M200 fall monitors, trusted by healthcare facilities for over 20 years. Sleek and durable, the multi-featured M210 offers a single, easy-to-configure platform to support your unique fall management program.

How It Works

With the M210 Fall Management System, you can connect two pressure-sensitive pads to a single monitor (one on a bed, another on a side chair or wheelchair). Select a time delay setting that best protects each person in your care. When a patient starts to exit the bed or chair, the monitor issues visual and auditory alarms to alert staff that immediate assistance is required. If a nurse call system is connected, alarms can be sent to the nurse call station.

The M210 monitor is easy to configure with a mobile app, allowing you to wirelessly configure the fall monitor to the individual needs of each patient or resident. Choose the configuration that best suits the criteria to successfully execute care plans and program goals: electric or battery power, auto-reset or non-auto-reset, audible or silent alerts, and tamper-resistant or easy-access rear cover.

The M210 Fall Monitoring System integrates with virtually any nurse call or e-call solution—to become an indispensable part of your fall management program.



Simplify Configuration and Deployment

The M210 Fall Monitor can be easily configured and firmware updated using the mobile app.

Long-Term Partnerships, Long-Term Value

Exceptional care requires teamwork, which is why we select only quality partners to deliver great products and service at the best price. Each of our trusted providers offer an unmatched level of personalized service and support to ensure end-to-end success. We value our partners and further stand by them with additional training to ensure excellent customer service and support at each stage of implementation.

STANLEY
Healthcare

About STANLEY Healthcare

Over 15,000 hospitals and senior living communities rely on STANLEY Healthcare solutions to empower caregivers to deliver better care. STANLEY Healthcare is a part of Stanley Black & Decker and a proud supporter of the Alzheimer's Association®. Learn more at stanleyhealthcare.com.

75 Portsmouth Blvd., Suite 220, Portsmouth, NH 03801 | 1-888-622-6992 | stanleyhealthcare.com | stanleyhealthcare@sbdinc.com